


Lead positions and one-channel rule used by the model

A. Realistic MLII torso electrode placement



B. What each model actually receives

NORMAL AUTOENCODER
18 NSRDB subjects
first stored channel only
WFDB name: ECG1
placement not documented

ABNORMALITY DETECTOR
MITDB MLII only
MLII selected by signal name
record 114: MLII is channel 2
records 102/104 excluded

SUBTYPE CLASSIFIER
same MLII window
same 242-value feature vector
same four autoencoder residuals
no second ECG lead is used

Therefore the performance numbers describe MLII temporal continuation, not a two-lead or unseen-patient result.