



# Real ECG cleaning for the two one-channel model inputs

**A. Autoencoder: NSRDB record 16265, first channel (ECG1)**


**1. Raw source signal**



**2. Resample to 128 Hz + 0.5--50 Hz band-pass**




**3. Standardise each non-overlapping 4 s window**




**B. Detector/classifier: MITDB record 100, MLI1**

**1. Raw source signal**



**2. Resample to 128 Hz + 0.5--50 Hz band-pass**



**3. Standardise each R-aligned 4 s beat window**

