


Chronological split within each MITDB record


Percentages refer to accepted-beat order, not elapsed recording time.

A Select model settings



Validation chooses the model candidates and the abnormality threshold.

B Refit, freeze, then evaluate



The model and threshold remain fixed while scoring the final test segment.

Hatched markers: 16 beats excluded immediately before each marked boundary. Markers are enlarged for visibility. The 64% gap is restored during final refitting.

Interpretation: later beats from known patients; the test had been seen in earlier project development, so this is a retrospective evaluation.